FR

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EN

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, the basic safety precautions below should always be followed.

READ ALL INSTRUCTIONS BEFORE USING.

DANGER:

To reduce the risk of electrocution:

- Do not handle plug with wet hands.
- Do not immerse in water or other liquid.
- Do not use while bathing or shower stall.
- Do not reach for a product that has fallen into any liquid. Unplug immediately.
- Check the power supply cord for damage before the first use and during the life of the product.

WARNING:

- To reduce the risk of burns, electric shock, fire, serious injury or damage to tissue:**
- Do not plug this device into a voltage system that is different from the voltage system specified on the device.
 - Do not use this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into any liquid.
 - Contact your local distributor or Broxo, which can be found at www.broxo.com.
 - Do not direct water under the tongue, into the ear, nose or other delicate areas. This product is capable of producing pressures that may cause serious damage in these areas. See operating instructions for correct usage.
 - Use this product only as indicated in these instructions or as recommended by your dental professional.
 - Only use tips and accessories recommended by Broxo.
 - Keep electrical cord away from heated surfaces.
 - Do not operate where oxygen or aerosol sprays are being used.
 - Do not use iodine, bleach, or water insoluble concentrated essential oils in this product.
 - Use of these can reduce product performance and will shorten the life of the product.
 - Fill reservoir with water or other dental professional recommended solution only.
 - Do not use while wearing any oral jewellery. Remove prior to use.
 - Do not use if you have an open wound on your tongue or in your mouth.
 - If your physician or cardiologist has advised you to receive antibiotic premedication before dental procedures, you should consult your dentist before using this instrument or any other oral hygiene aid.
 - For house hold use only. Do not use outdoors.
 - Closely instruct and supervise children and individuals with special needs in the proper use of this product.
 - Do not use this product for more than 4 minutes in each 10 minutes period. Longer usage can overheat the motor and cause product failure as well as voiding the warranty.
 - Children should be supervised to ensure that they do not play with the appliance.

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