DEUTSCH

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ENGLISH

- INSTRUCTIONS FOR WALL MOUNTING**
1. For better ease of use, it is recommended to place the holder on the wall to the left of the user or washbasin.
 2. Place the bracket against the wall in the desired position with the 2 fixing holes pointing upwards, making sure that the 4 rubber feet of the device mounted on its support do not come into contact with any area. These 4 feet are intended only for the storage of it, without use.
 3. Mark the position on the wall of the holder's three mounting holes with a pencil or sharp point.
 4. Drill three holes, approximately 3/16" in diameter, into the wall.
 5. Insert one of the supplied plastic plugs in each one of the drilled holes.
 6. Position the holder so that the mounting holes are located over the plastic plugs and secure the holder onto the wall by inserting a screw in each of the plastic plugs.
 7. Make sure the bracket is securely fastened, and slide the BROXO® Orajets™ unit onto the holder using the slots provided on the backside of the unit.

- HOW TO USE THE BROXO® OraJets™ FRACTIONNATED IRRIGATING SYSTEM**
1. Fit the Jet tip nozzle to the handle
 2. Fill the reservoir with lukewarm water and position it firmly in place on the BROXO® OraJets™.
 3. Turn the right-hand regulator knob to position 4.
 4. Place the Jet tip in the mouth, angled (45°) slightly towards the gum line and aimed so that water will penetrate the spaces between the teeth. The tip should not quite touch teeth and gums.
 5. Turn the left-hand motor switch to position 1 (on).
 6. Adjust the right-hand regulator knob to the most appropriate setting, from 1 to 7. Thoroughly irrigate your teeth and gums for at least one minute, directing the flow of water along the gum line. To prevent missing any area, irrigate following a definitive sequence.
 7. To momentarily stop irrigating, press on the knob labeled "STOP" situated on the handle.
 8. It is recommended to keep irrigating until the entire reservoir has been emptied.
 9. Turn the left-hand motor switch to position 0 (off).

- HOW TO USE THE BROXO® OraJets™ TONGUE CLEANING SYSTEM**
1. Fit the "Scrape & Flush™" nozzle to the handle
 2. Fill the reservoir with lukewarm water and position it firmly in place on the BROXO® OraJets™.
 3. Turn the right-hand power regulator knob to position 4.
 4. Turn the left-hand motor switch to position 1 (on).
 5. Adjust the right-hand regulator knob to the most appropriate setting, from 1 to 7. Thoroughly clean your tongue. To momentarily stop irrigating, press on the knob labeled "STOP" situated on the handle.
 6. It is recommended to keep irrigating until the entire reservoir has been emptied.
 7. Turn the left-hand motor switch to position 0 (off).

- HOW TO USE THE THE BROXO® OraJets™ NEBULIZING SYSTEM**
1. Fit the nebulizer nozzle to the handle
 2. Fill the reservoir with lukewarm water and position it firmly in place on the unit.
 3. Turn the right-hand regulator knob no higher than position 4. If the nebulizing cloud takes too long to form, repeatedly drive the right-hand regulator from position 1 to 7 before settling on position 4; or remove the nebulizer nozzle, switch the device on, re-attach the nebulizer.
 4. Turn the left-hand motor switch to position 0 (off).

- PARTIAL DISMANTLING OF THE PARTS**
- The BROXO® OraJets™ unit comes equipped with a special patented system that simplifies cleaning by allowing the dismantling of the complete hose/instruments handle and the mobile tray.
- Instructions to dismantle the cleanable parts of the unit:
1. Insert a thin flat tool under the rotating connector, where the hose enters the unit through the mobile tray.
 2. Use the tool as an arm of lever to push out the rotating connector.
 3. The complete hose/instrument handle and the mobile tray assembly should now be easily removable.
 4. To reassemble the unit, proceed in the opposite direction by placing the mobile tray back in its location and then firmly pushing the rotating connector into the hole until it snaps back into place.

- CLEANING OF THE APPLIANCE**
- Except for the nozzle handle, which may be rinsed regularly under running water, the main unit shall only be cleaned by means of a wet cloth which has been moistened with lukewarm water (if necessary mixed with cleaning products suitable for plastic materials).

- DANGER:** Do not immerse the unit in water or rinse it under running water.
- If the unit is used where the water contains a high percentage of limestone, the appliance shall be cleaned from time to time by running through its hydraulic system a reservoir of white vinegar. After this operation the appliance shall be rinsed by running at least one full reservoir of water. Nozzles can also be clogged by scale. Cleaning of the apertures shall be made by using a toothbrush under running water and this after the central piece of the nozzle has been pushed back by pressing on it. If still clogged, the nozzle shall be immersed in white vinegar or a chemical descaling solution suitable for plastic materials. Always rinse thoroughly under running water after cleaning.

- SERVICING**
- Any other servicing should be performed by an authorized representative.

- STORAGE OF THE DEVICE**
- When the device is not in use and has been removed from the wall where it shall normally always be secured, it shall be stored in a clean place and the cord shall be rolled and protected against any damage (heat, mechanical, chemical).

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